

READY•STEEL™

Perfectly Sterile. Perfectly Maillefer.

The stainless steel Maillefer instruments you trust are now available in pre-sterilized blister packs for single-patient use.



**WE
KNOW
ENDO.**

Introducing Maillefer **READY•STEEL**™

Ready When You Are

DENSTPLY Maillefer is an advocate of raising medical standards and is introducing the first step to a new complete sterile stainless steel solution: **READY•STEEL**™. **READY•STEEL** files are pre-sterilized and come with perforated blister packs for single-patient use, to bring you peace of mind during your endodontic procedures.



Reduced risk of file breakage*

Endodontic files exhibit signs of wear during normal use, which increases dramatically with multiple use. Single use reduces the risk of file breakage and increases patient safety.



Optimal cutting efficiency*

Research has shown that using files for multiple cases reduces their cutting efficiency, and leads to a poor performance. Cutting efficiency is achieved when a new instrument is used for each patient.




No risk of cross-contamination*

Research shows that certain bacteria and tissue remnants may remain on the instrument, post-sterilization. Single-patient use will eliminate such risk.

*Data available upon request.



Product	Details	Specs
FlexoFiles® Stainless Steel 	Uses: Exploration, canal patency, establishing glide path, cleaning and shaping, and apical gauging Features: Markings on end of handle, enhanced flexibility, safe-ended non-cutting Batt® tip, 25 and 31 mm have x-ray visible calibration in rings, and silicone stoppers	Lengths: 18, 21, 25, 31 mm** Sizes: 06, 08, 10, 12, 15, 17, 20, 22, 25, 27, 30, 32, 35, 37, 40, 45-80, 90-110, 120-140 Pack of: 6, assorted available Cross Section: ■ in sizes 6, 8, 10 ▲ in larger sizes
Senseus® FlexoFiles® Stainless Steel 	Uses: Exploration, canal patency, establishing glide path, cleaning and shaping, and apical gauging Features: Markings on end of ergonomic handle, flexible, improved tactile feedback, 25 and 31 mm have x-ray visible calibration in rings, and silicone stoppers	Lengths: 21, 25, 31 mm** Sizes: 6, 8, 10, 15-60 Pack of: 6, assorted available up to size 80 Cross Section: ■ in sizes 6, 8, 10 ▲ in larger sizes
K-File Stainless Steel 	Uses: Exploration, canal patency, establishing glide path, cleaning and shaping, and apical gauging Features: Markings on end of handle, safe-ended non-cutting tip, and silicone stoppers	Lengths: 21, 25, 31 mm** Sizes: 6, 8, 10, 15-80 Pack of: 6, Assorted Available Cross Section: ■
C+File Stainless Steel 	Uses: Accessing and cleaning calcified root canals Features: Ultra-stiff design, markings on end of handle, silicone stoppers and calibration rings	Lengths: 21, 25 mm Sizes: 6, 8, 10, 15 Pack of: 6 Cross Section: ■
Hedstrom File Stainless Steel 	Uses: Probing, permeabilization or extraction of debris Features: Markings on end of handle, safe-ended non-cutting tip, and silicone stoppers	Lengths: 21, 25, 31 mm** Sizes: 10, 15-80 Pack of: 6, assorted available Cross Section: ●
K-Reamer Stainless Steel 	Uses: Exploration, canal patency, establishing glide path, cleaning and shaping, and apical gauging Features: Markings on end of handle, safe-ended non-cutting tip, and silicone stoppers	Lengths: 21, 25, 31 mm** Sizes: 6, 8, 10, 15-60 Pack of: 6, assorted available up to size 80 Cross Section: ▲

**Not available in all sizes. Please contact your sales rep for details.

Why Pre-sterilization?

READY•STEEL™ files are conveniently packaged in pre-sterilized blister packs. In addition to its convenience, using a new pre-sterilized file for each patient will:

- Control costs by doing away with cleaning and sterilizing files.
- Maintain optimal cutting efficiency with a new file for each patient.
- Eliminate risk of cross-contamination from inconsistent sterilization procedures.



DENTSPLY MAILLEFER / T: 1-800-924-7393 / F: 1-800-924-7389 / MAILLEFER.COM